

10x and cellranger


COURSE ON SCRNA-SEQ DATA ANALYSIS

Fundamentals of 10x Chromium and preprocessing

Jan Kubovčíak · 7.-10. September 2026


with acknowledgments to Deepak Tanwar and SIB

Institute of Molecular Genetics of the Czech Academy of Sciences




How to capture cell of origin?

- Gel Beads-in-emulsion (GEM)
 - Way to label mRNA from the same cell





How to deliver the label?


- 10x chromium instrument




10,000+ GEMs in 4 minutes.



All captured transcripts from single cell:
identical + unique barcode



- reverse transcription
- breaking GEMs
- fragmentation
- primer ligation
- index PCR



Illumina sequencing

Sequencing output

- Filenames notation

sample id	lane	read type	
ETV6-RUNX1_1_S1	L001	R1	001.fastq.gz
ETV6-RUNX1_1_S1	L001	R2	001.fastq.gz
ETV6-RUNX1_1_S1	L001	I1	001.fastq.gz



- Files from inside



ETV6-RUNX1_1_S1_L001_R1_001.fastq.gz

@SRR9264343.187191690 K00271:135:HK3J7BBXX:3:2110:12591:47964 length=26

GGCAATTCATTGGGCCCCGAAGCTC

+

AAF-FJFJJJJJJJJJJJJJJJJJJ

@SRR9264343.5519124 K00271:135:HK3J7BBXX:3:1112:2950:8295 length=26

GGCAATTCATTGGGTGAGGTTGTAGG

+

AAFFJJJJJJFFJFJJJJJJJJJJ

@SRR9264343.149664941 K00271:135:HK3J7BBXX:3:2217:29579:19267 length=26

TCGGATGTCTGAAAGAGTGCTAGAAT

...

ETV6-RUNX1_1_S1_L001_R2_001.fastq.gz

@SRR9264343.187191690 K00271:135:HK3J7BBXX:3:2110:12591:47964 length=98

GTAGTAAAAGACTGGTTAATGATAACAATGCATCGTAAACCTTCAGAAGGAAAGGAGAATGTTTTGTGGACCACTTCGTTTTCTTTTTGCGTGTG

+

-777<<A-<FF<<A7-<FAFAFFJJJJJJFFFFJJFFJJFFJJFFJJJ<FJJAFJJFJJ<A-<AA<FFF)<A--7-7--AFFJ--AFafa7<AA-

@SRR9264343.5519124 K00271:135:HK3J7BBXX:3:1112:2950:8295 length=98

AGACTTTCTACCTGGTCATATACTCTGCAGCTGTTAGAATGTGCAAGCACTTGGGGACAGCATGAGCTTGCTGTTGTACACAGGGTATTTCTAGAAGC

+

--77-77--7-77<-7AFAJJJJJ7A7AFF7-A-7FFJJFA-AF77FF<-A<FFJJF-)FF7F<7FA7FAAAFJ-A--A<AAFJJ<F----7--<F77

@SRR9264343.149664941 K00271:135:HK3J7BBXX:3:2217:29579:19267 length=98

GGTCCGCAAGCAGGTGGTGCACATCCCGTCCTTCATTGTCCGCTGGCTTCCAGAAGCGCATCGACTTCTCTGCGCTCTCCCTACGGGGGTGGCC

...

ETV6-RUNX1_1_S1_L001_I1_001.fastq.gz

@SRR9264343.187191690 K00271:135:HK3J7BBXX:3:2110:12591:47964 length=8

CCTAGACC

+

AAFFJJJ

@SRR9264343.5519124 K00271:135:HK3J7BBXX:3:1112:2950:8295 length=8

CCTAGACC

+

AAFFFFJF

@SRR9264343.149664941 K00271:135:HK3J7BBXX:3:2217:29579:19267 length=8

CCTAGACC

...

Data preprocessing

1.Mapping

2.Quantification

3.Cell calling

- For 10x data use [cellranger](#)
- Alternatively [Alevin](#), [STARsolo](#)

Cellranger reference data

- Human, mouse and rat available for download
- Other organisms can be prepared with `cellranger mkref`
 - Use [Ensembl](#) for input data
 - Add exogenous sequences (eGFP etc.)
 - [Docs](#)

Counting UMIs vs alignments

- UMI: Unique Molecular Identifier:
 - Identifies each molecule (i.e. sequence) uniquely
- Molecules from a common PCR template -> carry the same UMI
- By counting UMI: correct for PCR replicates

Counting UMIs vs alignments




- Reads with **different UMIs** were derived from **different molecules** and are **biological** duplicates
- Reads with the **same UMI** originated from the **same molecule** and are **technical** duplicates – the UMIs should be collapsed to be counted as 1

	Cellular barcode	UMI		
Cell 1	TTGCCGTGGTGT	GGCGGGGA	CGGTGTTA	DDX51
	TTGCCGTGGTGT	TATGGAGG	CCAGCACC	NOP2
	TTGCCGTGGTGT	TCTCAAGT	AAAATGGC	ACTB
				
Cell 2	CGTTAGATGGCA	GGGCCGGG	CTCATAGT	LBR
	CGTTAGATGGCA	ACGTTATA	ACGCGTAC	ODF2
	CGTTAGATGGCA	TCGAGATT	AGCCCTTT	HIF1A
				
Cell 3	AAATTATGACGA	AGTTTGTA	GGGAATTA	ACTB ← 2
	AAATTATGACGA	AGTTTGTA	AGATGGGG	1
	AAATTATGACGA	TGTGCTTG	GACTGCAC	RPS15
				
Cell 4	GTTAAACGTACC	CTAGCTGT	GATTTTCT	GTPBP4
	GTTAAACGTACC	GCAGAAGT	GTTGGCGT	GAPDH
	GTTAAACGTACC	AAGGCTTG	CAAAGTTC	ARL1 ← 2
	GTTAAACGTACC	TTCCGGTC	TCCAGTCG	2
				
				
				

(Thousands of cells)

Cell calling

- Cell = high UMI content
- Empty droplet = low UMI content vs background



Cellranger count output

ETV6-RUNX1_1

Alerts

The analysis detected ⚠️ 1 warning.

Alert	Value	Detail
⚠️ Fraction of RNA read bases with Q-score >= 30 is low	59.4%	Fraction of RNA read bases with Q-score >= 30 should be above 65%. A lower fraction might indicate poor sequencing quality. This is Read 1 for the Single Cell 3' v1 chemistry and Single Cell 5' paired end, Read 2 for the Single Cell 3' v2/v3 chemistry and Single Cell 5' R2-only)

Summary [Analysis](#)

3,091

Estimated Number of Cells

68,259


Mean Reads per Cell

1,717

Median Genes per Cell

Sequencing ?	
Number of Reads	210,987,037
Number of Short Reads Skipped	0
Valid Barcodes	98.2%
Valid UMIs	100.0%
Sequencing Saturation	84.4%
Q30 Bases in Barcode	96.4%
Q30 Bases in RNA Read	59.4%
Q30 Bases in UMI	96.5%



Mapping ?	
Reads Mapped to Genome	95.8%
Reads Mapped Confidently to Genome	92.9%
Reads Mapped Confidently to Intergenic Regions	5.2%
Reads Mapped Confidently to Intronic Regions	25.5%
Reads Mapped Confidently to Exonic Regions	62.2%
Reads Mapped Confidently to Transcriptome	58.2%
Reads Mapped Antisense to Gene	1.2%



Sample	
Sample ID	ETV6-RUNX1_1
Sample Description	
Chemistry	Single Cell 3' v2
Include introns	False
Reference Path	...nger/refdata-cellranger-GRCh38-3.0.0
Transcriptome	GRCh38-3.0.0
Pipeline Version	cellranger-6.0.1

Cellranger count output

- Some tertiary analysis for EDA
- Advanced analysis with 10x Loupe Browser
- Or use third-party tools



Sequencing Saturation ?

Median Genes per Cell ?



10x portfolio for single cell



Flex

Fix, batch, and run on your schedule.

Protein coding gene coverage

Performs well with low-quality and FFPE samples

Multimic readouts from the same cell

Gene expression
Protein
CRISPR

Throughput options

Up to 8M cells (1–3,072 samples) per run

Products

[Flex >](#)



Universal

Species agnostic. Maximum versatility.

Whole transcriptome coverage

Delivers broadest set of information, including isoforms, SNPs, etc.

Multimic readouts from the same cell

3' or 5' gene expression
TCR/BCR
Protein
CRISPR

Throughput options

Up to 160,000 cells (1–8 samples) per run

Products

[Universal 3' >](#)
[Universal 5' >](#)



Epi Chromatin

Unmask epigenomic profiles.

ATAC-seq chemistry

Explore open chromatin regions
Link directly to 3' gene expression (Multiome kit)

Multimic readouts from the same cell

Chromatin accessibility
3' gene expression

Throughput options

Up to 80,000 nuclei (1–8 samples) per run

Products

[Epi ATAC >](#)
[Epi Multiome >](#)

Thank you for your attention!